Ordering Information

Part Number	Description
Intelect Neo Vacuum Module	
2765CS	Intelect Advanced Color Stim + EMG
2755CS	Intelect Color Stim
2773MS	Intelect Advanced Monochromatic Stim
Standard Accessories	
2771	sEMG and sEMG + Electrical Stimulation Module (REF 2765CS only)
27465	Patient Data Card (1x)
10648	Nylatex Wrap
79967	6 x 8 cm carbon electrodes (4x)
79970	6 x 8 cm sponges (4x)
42192	Dura-Stick® Plus 5 cm round electrodes (4x)
4248	Ultrasound gel (Combo only)
27455	User Manual on CD
Optional Accessories	
27469	Patient Interrupt Switch
2771	sEMG and sEMG + Electrical Stimulation Module [incl. Durastick 3 cm round electrodes (4x)]
2767	Battery Module
2766	Laser Therapy Module
2770	Channel 3/4 Stimulation Module
2785	Vacuum Electrode Module
2774	Vacuum Module w/Therapy System cart
27133	Intelect Advanced Carrying Bag
27508	Operator Remote (Ch 1/2)
27079	Operator Remote (Ch 3/4)
2775	Therapy System Cart
79977	HiVolt Probe
57007	Microcurrent Probe
114.610 EXP	Intravaginal Probe

For Electrodes, please see page 32

Technical Specifications

Mains Power:	100-240 VAC, 50-60 Hz
Weight:	Module only: 2.3 kg
Dimensions*:	Module only: 26 x 14 x 28 cm
Electrical Safety Class:	Class I, Type BF
Safety Tests:	IEC 60601-1, IEC 60601-1-2, IEC 60601-2-10

* Dimensions are listed as L x W x H

Intelect® Advanced Stim

chattanooga



Features

Two channel electrotherapy

The unique modularity found in the Intelect Advanced Therapy System allows for easy customization. Additional components that can be added to the Stim base now or later with a simple upgrade include: EMG, 2 additional channels of stimulation, laser therapy, vacuum, battery.

- A complete range of clinical waveforms
- 2 independent electrotherapy channels
- Over 200 Clinical Protocols
- Over 100 User-defined protocol slots
- Over 10 Quick-link Indications
- Clinical library includes anatomical and pathological libraries on vibrant LCD screen
- Store treatment sessions easily on Patient Data Cards equipped with the system
- Sequencing to direct a designated therapeutic course without having to modify unit parameters
- Constant current/constant voltage modes
- On screen pad contact quality indicator
- 3 year warranty (probes and accessories 1 year)

