LANDICE

SPORTS/PHYSICAL THERAPY & REHABILITATION

REHABILITATION TREADMILLS - ELLIPTICALS - BIKES

The Landice line of rehabilitation products provides the very best in functional usage and reliability. From reversing treadmills to low-impact ellipticals and easy-to-start self-powered bikes, Landice fitness equipment provides the perfect tools to assist a wide variety of patients in their recovery and rehabilitation.



REHABILITATION EQUIPMENT

L7 / L8 Rehab Treadmills
E9 Ellipticals
R9 / U9 Recumbent & Upright Bikes

THE ROAD TO RECOVERY STARTS HERE

L7 & L8 REHAB TREADMILLS

Made with the finest materials, the L7 and L8 Rehabilitation Treadmills provide the best and safest support for all your rehabilitative programs. With a low step-up height, zero starting speed, 300 micro-amp isolation leakage kit and extended parallel medical handrails, these treadmills are ideal for any physical therapy and sports conditioning center. A number of unique options are also available, such as a reversing treadbelt for downhill walking, 100 micro-amp hospital isolation leakage kit, and an orthopedic suspension system. For a complete list of options please visit us at Landice.com.

TREADMILL FRAME COMPARISON

Solid aluminum construction - rust-free and durable

	L8	L7
User Capacity	500 lbs.	400 lbs.
Equipment Weight	305 lbs.	265 lbs.
Dimensions	35″ W x 83″ L x 61″ H	35″ W x 76″ L x 59″ H
Running Surface	22″ W x 63″ L	20″ W x 58″ L
Power	110 VAC, 60Hz, 15A 220 VAC, 60Hz, 15A (optional)	110 VAC, 60Hz, 15A 220 VAC, 60Hz, 15A (optional)
Rollers	3 1/2″ 22 lbs.	2 1/2″ 14 lbs.
Step-up Height	7 1/2″	5 1/2″
Speed	0.1 to 12 mph	0.1 to 12 mph
Incline	0% to 15%	0% to 15%
Drive System	4 HP Continuous Duty Motor	4 HP Continuous Duty Motor
Deck	1″ - Reversible	1″ - Reversible
Control Panel	Rehabilitation	Rehabilitation
Warranty*	5 yrs. parts/2 yrs. wear items/ 1 yr. labor	5 yrs. parts/2 yrs. wear items/ 1 yr. labor

*Treadmills are covered by the LTD warranty which offers coverage for usage of less than 5 hours per day.



SUPERIOR REHAB EQUIPMENT

COMFORTABLE AND SAFE REHAB

BIKE SPECIFICATIONS

	R9	U9
Frame Type	Step-through recumbent with low step-over height	Upright
Pedal Type	Self-balanced pedals with adjustable straps	Self-balanced pedals with adjustable straps
Seat Adjustment	Forward/backwards with 15 settings	Up/down with 12 settings
Flywheel Weight	15 lbs.	15 lbs.
Resistance System	Eddy current magnetic braking system	Eddy current magnetic braking system
Resistance Range	1-20	1-20
Handlebar	Console grip bars, seat grips with contact heart-rate monitoring and resistance controls	Racing style handlebars with contact heart-rate and resistance controls
Display	LED with multiple feedback windows	LED with multiple feedback windows
Heart Rate	Contact heart rate for heart rate control/ wireless chest-strap (optional)	Contact heart rate for heart rate control/ wireless chest-strap (optional)
Accessories	Water bottle holder, reading rack, accessory tray	Water bottle holder, reading rack, accessory tray
Cooling Fan	3 speeds	3 speeds
Power Requirement	Self-powered	Self-powered
Transport Wheels	Yes	Yes
Levelers	Yes	Yes
Assembled Weight	205 lbs.	159 lbs.
Dimensions	65.74" L x 25.31" W x 52.20" H	44.56" L x 24.72" W x 58.11" H
Max User Weight	350 lbs.	350 lbs.
Warranty*	5 yrs. parts/1 yr. labor/90 days for high-wear items	5 yrs. parts/1 yr. labor/90 days for high-wear items

*Bikes are covered by the Club warranty which offers coverage for usage of more than 5-hours per day.

R9 / U9 BIKES

The Landice R9 & U9 bikes offer a comfortable and safe rehabilitation option. The wide step-through design on the R9 allows for easy access along with superior back support. The U9 offers well-cushioned seating with adjustable height levels to support a personalized lower-body rehabilitation regimen.



NATURAL BIOMECHANICS - LOW IMPACT

E9 ELLIPTICAL TRAINER

Patients at any athletic level from beginner to professional will find physical benefits and rehabilitative functionality on a Landice.

The E9 offers one of the only fitness units that enables both upper and lower body exercise capabilities. Patients can stand upright on the frame's stationary footpads and perform an isolated workout to strengthen the upper body. This back-and-forth motion is considered a preferable position as opposed to a seated rotational movement.

For the lower body, the center-drive technology replicates the biomechanics of walking or running. The E9 provides a low impact alternative to treadmills.

The elliptical also features an accessible, safe, stable and patient-friendly exercise platform. Users can easily grab the side rails, step-up and safely start a workout. The pedals will move only when the patient engages the upper body handrails and/or foot pedals.



E9 ELLIPTICAL FEATURES

User Capacity	500 lbs.	**Stride Length**	Adjustable 17" - 23"
Active Footprint	35" W x 76" L x 67" H	**Drive System**	Center Drive
Physical Footprint	35" W x 68" L x 67" H	**Flywheel**	11" diameter, 20 lbs.
Power	110 VAC, 60Hz, 15A 220 VAC, 60Hz, 15A (optional)	**Panels**	Pro, Pro Sports, Cardio, Executive
Handrails	2" Medical-Grade Side Rails	**Pedals**	2" Spacing with Orthopedic Gel Inserts
Warranty*	5 years parts/2 years wear items/1 year labor		

*Ellipticals are covered by the Club warranty which offers coverage for usage of more than 5-hours per day.

The Landice E9 is available in 4 control panel styles: Executive, Cardio, Pro Sports, and Pro. For more information on these consoles, please visit us at Landice.com.

Landice engineers, assembles and tests the finest specialty fitness equipment on the market and is one of the only manufacturers to do so in the USA. Landice has been recognized repeatedly by leading consumer rating publications as having the world's finest treadmills and ellipticals.

Visit us at **www.landice.com** for additional options and complete warranty information.