Laser Module



Features

Low Level Laser Therapy

- Fully functional device with independent control over all parameters
- Display of energy in choice of joules or joules/cm²
- Continuous and pulsed treatment modes with fully variable pulsed frequencies
- Real time feedback of energy delivered
- Full selection of clinical indications specific to laser applicator selected
- Wide variety of LED, Laser Diodes, & Cluster Probes available
- User-defined protocol slots on board
- Easy connection to existing Intelect Advanced Electrotherapy or Combination device
- Built in acupuncture point finder
- 3 year warranty (probes and accessories 1 year)

⚠ Visible and invisible laser radiation.
Avoid direct exposure to beam.
Class 3B Laser Product.



Technical Specifications	
Weight:	1 kg
Dimensions*:	26 x 21 x 4 cm
Electrical Safety Class:	Class 1, Type B
Safety Tests:	IEC 60601-1, IEC 60825-1, IEC 60601-2-22

* Dimensions are listed as L x W x H

Ordering Information	
Part Number	**Description**
Laser Module	
2766	Laser Therapy Module
Standard Accessories	
27525	Laser Protection Glasses (2x)
27455	User Manual on CD
Optional Accessories	
27900K	Remote Interlock Kit
Laser Applicators	
Singles	
27799	670nm LED Diode 10mW
27803	850nm Laser Diode 40mW
27840	850nm Laser Diode 100mW
27804	850nm Laser Diode 150mW
27841	850nm Laser Diode 200mW
27805	820nm Laser Diode 300mW
9 Diode Cluster	
27810	9 Diode Cluster Applicator, 290mW Total: 5 x 850nm 50mW Lasers, 4 x 670nm 10mW LED
27811	9 Diode Cluster Applicator, 540mW Total: 5 x 850nm 100mW Lasers, 4 x 670nm 10mW LED
27812	9 Diode Cluster Applicator, 1040mW Total: 5 x 850nm 200mW Lasers, 4 x 670nm 10mW LED
13 Diode Cluster	
27813	13 Diode Cluster Applicator, 265mW Total: 3 x 850nm 50mW Lasers, 7 x 670nm 10mW LED, 3 x 950nm 15mW LED
27814	13 Diode Cluster Applicator, 415mW Total: 3 x 850nm 100mW Lasers, 7 x 670nm 10mW LED, 3 x 950nm 15mW LED
27816	13 Diode Cluster Applicator, 715mW Total: 3 x 850nm 200mW Lasers, 7 x 670nm 10mW LED, 3 x 950nm 15mW LED
19 Diode Cluster	
27815	19 Diode Cluster Applicator, 325mW Total: LED only, 6 x 670nm 10mW LED, 7 x 880nm 25mW LED, 6 x 950nm 15mW LED
33 Diode Cluster	
27809	33 Diode Cluster Applicator, 565mW Total: LED only, 12 x 670nm 10mW LED, 13 x 880nm 25mW LED, 8 x 950nm 15mW LED
27802	33 Diode Cluster Applicator, 690mW Total: 5 x 850nm 50mW Lasers, 12 x 670nm 10mW LED, 8 x 880nm 25mW LED, 8 x 950nm 15mW LED
27807	33 Diode Cluster Applicator, 940mW Total: 5 x 850nm 100mW Lasers, 12 x 670nm 10mW LED, 8 x 880nm 25mW LED, 8 x 950nm 15mW LED
27808	33 Diode Cluster Applicator, 1440mW Total: 5 x 850nm 200mW Lasers, 12 x 670nm 10mW LED, 8 x 880nm 25mW LED, 8 x 950nm 15mW LED