Vacuum

chattanooga



Features

Vacuum Electrode Module

- 2 channel module (4 electrode connections) that can support 4 channels with the addition of a second module in the Therapy System Cart
- Quick control knobs make adjustment of cup strength simple
- Unique positioning of vacuum in Therapy Cart allows for efficient draining into resevoir below and makes reservoir bin cleaning easily accessible
- Easy connection to any Intelect Advanced Therapy System or Intelect Mobile Device
- Various sizes of Vacuum Electrode Cups available
- 3 year warranty (accessories 1 year)



Ordering Information

Part Number	Description
Vacuum Electrode Modules	
2774	Vacuum Electrode Module w/ Therapy System Cart
2785	Vacuum Electrode Module (not shown)
Standard Accessories	
56852	60 mm Vacuum Electrodes Stainless Steel (4x)
56806	60 mm Vacuum Sponges (4x)
56902	Ch 1/3 Vacuum Lead Hose (Black)
56901	Ch 2/4 Vacuum Lead Hose (Gray)
27455	User Manual on CD
Optional Accessories	
27657	30 mm Vacuum Electrodes Stainless Steel (4x)
27658	30 mm Vacuum Sponges (4x)
56899	60 mm Vacuum Electrodes Aluminum (4x)
56851	90 mm Vacuum Electrodes Stainless Steel (4x)
56900	90 mm Vacuum Electrodes Aluminum (4x)
56807	90 mm Vacuum Sponges (4x)
56902	Channel 1/3 Vacuum Lead Hose (2x) (Black)
56901	Channel 2/4 Vacuum Lead Hose (2x) (Gray)

Technical Specifications

Weight:	2 kg
Dimensions*:	23 x 13.7 x 24 cm
Electrical Safety Class:	Class I, Type BF
Safety Tests:	IEC 60601-1, IEC 60601-1-2

* Dimensions are listed as L x W x H